

## Specification

Microscope Head	Binocular, 30° inclined, 360° rotating, with tablet
Microscope Eyepieces	WF 10x/20
Microscope Nosepiece	Quadruple, reversed
Microscope Objectives	N-PLAN 4x, 10x, 40x, 400x (oil)
Microscope Stage	Double layer, 150x139 mm, moving range 75x33 mm
Digital Camera Resolution	3 MP
Operating System	Windows 10 64-bit
Image Capturing Software	OPTIKA Vision lite
CPU	Intel® Atom™ Z830, Quad core
LCD Display	LED 10.1" IPS Multi Touch Screen
LCD Resolution	1920x1080, 16/9 Full HD
Storage	Hdd 64 GB
Network	Wireless, Bluetooth 4.0
Control Buttons	Auto rotate off, volume control
Battery Technology	Lithium-ion battery
Also Included	Instructional manual, keyboard with touchpad, touch pen